Deutsch

English

Français

Gebrauchsanweisung
Operating Instructions
Mode d'emploi

Variotherm plus

CE 0124

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Table of contents

Page

1.0	Introduction	
1.1	Notes on operating instructions.....	5
1.2	Intended use	7
1.3	Function	7
1.4	Explanation of symbols	9
2.0	Safety advice	12
2.1	Advice for the user for maintaining hygienic status.....	13
3.0	Installation and start-up	17
3.1	Standard delivery	17
3.2	Illustrations	18 - 21
3.3	Connections	23
3.3.1	Electrical connection	23
3.3.2	Connecting a nystagmograph	23
3.3.3	Connecting a potential balance.....	23
3.3.4	Water connection	23
3.3.5	Connecting the irrigation handle	25
3.4	Start-up	25
4.0	Operation	
4.1	Adjusting temperatures	27
4.2	Selecting temperature levels.....	27
4.3	Adjusting stimulation time	27
4.4	Adjusting flow quantity	27
4.5	Description of operating modes	29
4.5.1	Thermal process for reducing bacterial count.....	29
4.5.2	Irrigation mode	29
4.5.3	Stimulation mode	31
4.5.4	Energy saving mode	31
5.0	Cleaning	
5.1	General information on cleaning and disinfection	33
5.2	Recommended disinfectants.....	35
5.3	Cleaning method for handle with jet connection for water irrigation tip.....	38
6.0	Maintenance and service	
6.1	Functional test.....	43
6.2	Visual inspection	44
6.3	Electrical test.....	44
6.4	Decalcification	44
6.5	Replacement of filter	45
7.0	Trouble shooting	50
8.0	Accessories and spare parts	
8.1	Accessories.....	53
8.2	Spare parts.....	53
9.0	Specifications	56
10.0	Disposal	58
11.0	Notes on EMC	62-64

General Standard Terms and Conditions

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1.0 Introduction

1.1 Notes on operating instructions

- The operating instructions contain important information on how to operate the Variotherm plus safely, correctly and effectively. The operating instructions are intended for new operating personnel and also should be continued to be used as a reference manual. Reading and understanding the operating instructions will assist the user in avoiding risks, reducing repair costs, lowering down-time, increasing reliability and extending the service-life of the equipment. **The operating instructions must always be kept available and near the equipment.** Prior to first use, please read the chapter "Safety advice", to be prepared for any possible dangerous situations.

The basic principles are:

Judicious and careful work provides the best protection against accidents!

Operational safety and readiness for use depend not only on your capabilities, but also on **care and maintenance** given to the Variotherm plus. Regular cleaning and service work are a must. Major maintenance and repair work should be carried out only by expert personnel authorised by ATMOS. For repairs, only use ATMOS spare parts to preserve the equipment warranty and ensure operational safety, readiness for work and to maintain the value of the equipment.

- The product Variotherm plus bears CE Marking CE-0124 according to the EEC guideline of the council for medical products 93/42/EEC and meets the basic requirements of this guideline.
- The quality management system applied at ATMOS has been certified according to international standards EN ISO 9001 and EN 13485.
- Reprints, and text extracts are allowed only with written permission from ATMOS.

Abbreviations / symbols in these operating instructions:

- Indicating a list
 - Subdivision of a list/activity

The recommended sequence must be followed in each case!

 Indicating particularly important advice!

 Describing the effect of an activity.

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1.2 Intended use

Main function:	A user-friendly compact device for the irrigation of the auditory canal and the stimulation of the vestibular organ
Application:	For the application on humans
Specifications of the main functions:	<p>Produces a jet of water for the ear irrigation at body temperature (37°C) with a flow rate of 400 ml/min for the removal of cerumen</p> <p>Stimulation of the vestibular organ with a defined flow and a fixed temperature for the functional testing</p>
Application organ:	Auditory canal to the drumhead
Duration of application:	For the short-term application on patients
Application environment:	Application environments are clinics and practices of ENT physicians and phoniaticians. The therapy with the ATMOS irrigation devices may only be performed by trained staff.
Contraindications:	Do not apply to an infected resp. contaminated auditory canal.

1.3 Function

- Operating the main switch on the examination unit automatically activates the thermal process for reducing bacterial count (see section 4.4.1) .
- Then, automatic change into the energy saving mode.
- When removing the irrigation handle from its holder, automatic change into the irrigation mode. Here, the auditory canal can be rinsed with water warmed up to 37°C and a flow of at least 400 ml/min.
- Possibility to change into the stimulation mode for stimulating the vestibule with a reduced flow of water. The **Variotherm plus** is equipped with a timer for preselecting the stimulation time.

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1.4 Explanation of symbols

	Pay attention to operating instructions ! as to ISO /7000/0434 DIN 30600/1008 IEC 348
	Type B equipment as to IEC 417
	Fuse as to IEC 417/5016, DIN 30600/0186
°C	Temperature in degree centigrade
S	Timer adjustment in seconds
	Start
	Stop
	Timer
	Cold stimulation level
	Warm stimulation level
	Irrigation level (water at a temperature value of 37°C)
	Heating ON
	Heating OFF (energy saving mode)
	Control output for connecting a nystagmograph (graphical recorder as to DIN 30600, IEC 417 5192)
	Equipotentiality DIN 30600 495, ISO 417 5021
	Connection for waste water
	Connection for water
Low	Reduced flow for stimulation of the vestibule
High	High flow for rinsing the auditory canal

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2.0 Safety advice

- The Variotherm plus unit is produced according to IEC 601/EN 60601 and listed in the following classes:
 - VDE Class of protection 1
 - Class IIa (EEC 93/42).
- The equipment must only be connected to a power outlet with a correctly installed ground/earthing contact according to local and state regulations.
- The unit must be installed by a specialist authorized by ATMOS (see chapter 3.3).
- The Variotherm plus may be used in **supervised operation** by qualified personnel which have been authorised by ATMOS and which have been trained in the operation of the equipment (IEC 601-1/EN 60601-1).
- The mains voltage specified on the type plate must match the power supply system.
- Prior to every use of the equipment, verify that the device is technically safe and in proper working condition. **Damaged leads** must be replaced immediately!
- Correct configuration in assembly of country-specific connections

green/yellow:	protective conductor (PE)
blue:	neutral conductor (N)
black resp. brown:	phase (L)
- It is essential that the country-specific requirements for connection of medical equipment to the public drinking water supply are considered.
When in doubt, please contact your local ATMOS partner.
- The control panel must be visible and in reach of the operator. Ensure the equipment is stable prior to use.
- ☞ ● In the thermal process for reducing bacterial count hot water is being ducted through the irrigation handle. Please do **not** take the handle from its holder or spray water!
- ☞ ● Metal parts can be hot!
- ☞ ● Prior to spraying, the water temperature must be checked by the user!
- Switch off main power switch after use and **close the water valve** of water supply.
- The Variotherm plus should only be operated in rooms used for medical purposes, and never in areas subject to explosion hazards and in oxygen rich environments.
- All additional equipment, which is connected to the analog and digital interfaces of the unit, must meet the requirements of relevant EN specifications (for inst. EN 60950 for data processing equipment and EN 60601 for electrical medical appliances). In addition, configurations must satisfy system specification EN 60601-1-1. When additional equipment is connected to the signal input or signal output section on the unit, the person making the connection is deemed "a system configuration operator" and as such is responsible for meeting the requirements of system specification EN 60601-1-1. For answers to additional questions, please contact your local specialist supplier or ATMOS technical service.
- The irrigation jet must not get into contact with contaminated material.
- Only use the warm water jet when a tube tip is installed!
- Take caution to avoid injury to the eardrum when introducing the water jet !
- For hygienic reasons jet connections must be changed after each patient. This also prevents retrograde germination of the warm water system.
- Use only for irrigation of auditory canal !
- ATMOS is not liable for personal injury and damage to property if
 - original ATMOS parts are not used,
 - the operating instructions are not completely followed,
 - assembly, new settings, alterations, extensions and repairs have been carried out by personnel not authorised by ATMOS.
- Please note:
A medical insulating transformer with earth leakage monitor or any similar safety system acc. to EN 60 601-1 is required, if several devices are connected over one common power supply. The transformer must correspond to the power consumption of all the devices to be connected.
- The flow rate has to be checked after switch on or at least once a day with a measuring cup. The quantity of 500 ml/min may not be exceeded. The water jet has to be straight.

2.1 Important instructions for the maintenance of the hygiene status of warm water irrigation units

For loosening cerumen in the auditory meatus and for the stimulation of the labyrinthine system, ATMOS offers the warm water units Hygrotherm plus (37°C) and Variotherm plus (20° - 47°C).

These units heat the drinking water which comes from the household connection to the preselected temperature. The operating instructions must be followed, or the increased number of pathogens in the rinse water can affect the health of sensitive patients.

Requirement for taking into use

- The water which is provided by the household connection must at least meet the WHO guidelines or the country-specific guidelines for drinking water.

Connection

- Fresh water must flow through the wall connection for a period of approx. 1 minute prior to the unit being connected.*
- Ensure that the installation is hygienically carried out! Disinfect the threaded connections with 70 % alcohol* before screwing them into place.
- Wear disposable gloves* when carrying out maintenance work on parts which come into contact with the water.
- Carry out thermal cleaning at least 3 times before patient use.*

* This work is carried out by the ATMOS service technician.

Running operation

- Before beginning the surgical operation, switch on the unit and await the thermal cleaning process.
- After the device had been idle for a long period of time, (weekend, during holidays, etc.) the thermal cleaning process must be initiated and repeated.
- We recommend that the device be switched off and on every 2 hours and not longer than 4 hours to start the thermal cleaning mode.
- Exchange and disinfection required for contact parts such as the jet connection after each patient (prevention of retrograde germification).
- Before every use, spray water and check the temperature.

Service

- Ensure that the suspended particle filter is hygienically changed (see commissioning) (dispose of filter and carefully clean the filter glass). Ensure that no contact is made with the contaminated parts.
- Carry out thermal cleaning after each service.
- We recommend that the number of germs in the water be regularly determined at the input and output of the unit. In the case of a complaint, data measurements must be prove whether the contaminants are "water" or "skin" bacteria.

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3.1 Standard delivery

Quantity

1	Variotherm plus basic unit
1	Irrigation handle
1	Double-barrelled hose
3	Jet connection short, straight (80 mm)
1	Splash protection
1	Rubber bush for splash protection
2	Three-holes gasket
1	Hose tips for jet connection (30 Stk)
1	G3/4i-G1/4a-screwing (water-tap G3/4 on filter G1/4)
1	Cap nut G3/4
1	Water filter complete
2	Gasket (f.G1/4a) 13x18x1
1	G1/4a-G3/4a-reducer (filter G1/4 on water hose G3/4)
1	Water supply hose G3/4i, L = 3 m
1	Waste water hose G1/4i, L = 3 m
2	Gasket (f. G1/4i) 10x15x1
1	Cardboard 455 x 340 x 250 mm
1	Padding for cardboard
1	Operating instructions

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D

Fig. 1. Front view

- ① Main switch
② Key switch f. heating ON/OFF (energy saving mode)
③ Coding switch for warm-stimulation level
④ Coding switch for cold-stimulation level
⑤ Coding switch for stimulation time
⑥ Key switch for selecting the irrigation level (37°C)
⑦ Key switch for selecting the warm-stimulation level (e.g. 44°C)
⑧ Key switch for selecting the cold-stimulation level (e.g. 30°C)
⑨ Key switch for start/stop of the stimulation
⑩ Temperature display (two-figure number, resolution 1°C), current value indication
⑪ Display of stimulation time (two-figure number, resolution 1s)

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D

Fig. 2. Rear view

- ① Type plate
- ② Control output for a nystagmograph
- ③ Unit plug with fuse compartment
- ④ Equipotentiality connection

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D

Fig. 3. Side view


- ① Support for irrigation handle
- ② Connection for double-barrelled hose
- ③ Connection for waste water hose
- ④ Connection for water supply hose

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D

Fig. 4. Irrigation handle

- ① Handle
- ② Connection for double-barrelled hose
- ③ Three-holes gasket
- ④ Adjusting screw
- ⑤ Set screw
- ⑥ Valve plunger
- ⑦ Jet connection
- ⑧ Splash protection
- ⑨ Cap
- ⑩ Double-barrelled hose

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The sprayer tube must be exchanged after each patient.

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UK

3.3.1 Electrical connection

- According to the directions of VDE 0107 and VDE 0100, medically used rooms have to be equipped with a *leakage current protective circuit* (FI protective circuit) with a nominal leakage current of < 0.03 A. Installation must correspond with VDE 0107.
- Connect the power cable to the socket on the unit (③, fig. 2).
- Insert the power plug in a correctly installed socket with a ground/earthing contact.

3.3.2 Connecting a nystagmograph

- When controlling an ENG (electro-nystagmograph) or a CNG (computer-nystagmograph) at output (②, fig. 2), please connect only recording equipment approved by ATMOS. Connecting cable available from ATMOS (see section 8.0).
- At the end of the stimulation time, a trigger signal for a nystagmograph is available at the rear 3-pole DIN socket. This outlet is an electrically isolated photo transistor. The emitter connection which has to be connected to the ground line of the trigger inlet at the nystagmograph is attached to pin 3. The open-collector connection is located on pin 1 and is to be connected to the positive supply voltage (+5 V or +12 V) by means of a pull-up resistance. The maximum collector voltage of the transistor may not exceed 80 MA.

3.3.3 Equipotentiality connection

- Connection for potential compensation (④, fig. 2). Connecting cable available from ATMOS (see section 8.0).

3.3.4 Water connection

- Local requirements:
 - Water-tap with G3/4" external thread.
 - Drinking-water!
 - Input water pressure: +2000...+6000 hPa.
 - Feed temperature: +7...+19.5°C (recommended), but at least 0.5°C below the required lowest cold stimulation value.
 - Drainage fittings (Art.No. 502.0880.0)
- Connection to water supply and waste water:
 - The country-specific connection requirements for the connection to the public water network have to be observed. To comply with the rules according to EN 1717, an unpressurized water separation with overflow is required. Therefore, ATMOS recommends to install the water separation unit ATMOS Aqua Clean, REF 502.1200.0 between the water supply and the device in order to meet the connection requirements according to EN 1717.
 - Prior to connecting the unit to the water supply the feed line must be flushed clean by opening the water tap for a minute and allowing free flow of the tap water.
 - Connect the filter unit to the water-tap by means of adaptor and gasket.
 - Insert the gaskets in the cap nuts of the water supply hose.
 - Join the nuts with the filter connection and the unit connection (④, fig. 3).
 - Insert the gaskets in the cap nuts of the waste water hose
 - Join the nuts with the unit connection (⑤, fig. 3) and the drainage fittings.
 - The water system is not equipped with any calcification protection. For drinking water with a water hardness of 14-21°d resp. 2.5-3.8 mmol/l (= hard water) and up from 21°d resp. up from 3.8 mmol/l (= very hard water) a calcification protection system is required. Please contact your local water supplier and/or plumber.

Water hardness	Millimol per litre	°dH	
1 (soft)	≤1.3	≤7.3	Calcification protection system is not required
2 (medium)	1.3 - 2.5	7.3 - 14	
3 (hard)	2.5 - 3.8	14 - 21.3	Calcification protection system is required
4 (very hard)	>3.8	>21.3	

Close the water tap when the unit is not in use!

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3.3.5 Connecting the irrigation handle

- Insert a "three-holes gasket" (③, fig. 4) in the double-barrelled hose and screw it (⑩, fig. 4) on the unit connection (②, fig. 3).

☞ Use only "three-holes gaskets" or the unit will not function correctly!

- Insert a "three-holes gasket" in the irrigation handle and screw it on the free end of the double-barrelled hose.
- Open the water-tap.
Check whether all connections are tight.

3.4 Start-up

- Insert the handle in its holder; the jet must point to the reverse side of the unit.
- Switch on the unit (①, fig. 1).
- The unit will display a test mode with digital numbers "8 8" and an audible signal.
- Automatic activation of the thermal process will begin for reducing bacterial count.
Duration: 5 min.
- When the thermal process is complete, the unit will switch to the irrigation mode (37°C).
- When the irrigation mode is not used (handle is not removed from its support) the unit change into the energy saving mode after 5 min.

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4.0 Operation

UK

- Before switching on the unit, make sure that the handle is in its holder (jet must point to reverse side of the unit).

4.1 Adjusting temperatures

- Number of temperature levels: 3
 - One level fixed to 37°C irrigation temperature.
 - Two variable temperature levels (20°C - 47°C) (48°C + 49°C only for testing purposes).
- Temperate setting by coding switches (3, 4, fig. 1).
 - left switch: for adjusting the "ten" partition
 - right switch: for adjusting the "one" partition
 - ↳ lower keys (+): Temperature increase
 - ↳ upper keys (-): Temperature decrease
- Standard settings:
 - Level for irrigation mode: 37°C fixed
 - Level for cold stimulation: 30°C
 - Level for warm stimulation: 44°C

4.2 Selecting temperature levels

- After operating the main power switch, the automatic activation of the thermal process begins for reducing bacterial count, followed by the automatic activation of the temperature level "irrigation mode" (37°C).
- For selecting the desired temperature level, use the respective key (6, 7, 8, fig. 1).
 - ↳ Display of the active level by illuminated diodes.
 - ↳ Display of the water temperature (current value) in °C.

4.3 Adjusting stimulation time

- By means of coding switch (5, fig. 1).

4.4 Adjusting flow quantity

- Prior to stimulation turn the set screw (1, fig. 6) upwards and arrest it.
- 2 adjusting possibilities:

High:	Flow quantity for rinsing the auditory canal
Low:	Flow quantity for stimulating the vestibular organ

Prior to ear irrigation the set screw is turned downwards, inserted into the handle and then adjusted by turning it downwards again.

- For fine adjustment of the flow the adjusting screw (2, fig. 6) is turned in respectively out.
 - ↳ Periodically check the water flow for the stimulation of the vestibular organ and readjust, if necessary.

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4.5 Description of modes

4.5.1 Thermal process for reducing bacterial count

Purpose:

Reduction of bacteria count in the warm water circuit to prevent or clear a contaminated system.

Activation:

Occurs each time the unit is switched on.

- Condition for starting the correct thermal process for reducing bacterial count:
 - The irrigation handle must be inserted in the holder on the side (jet pointing to reverse side of the unit).
 - ↳ If not, a warning signal is given (warning also sounds when the handle is removed during the thermal cleaning procedure).
 - ↳ Additionally, the error code "F0" is shown (disappears when handle is returned).
- ☞ • In the thermal process for reducing bacterial count, hot water is passing through the irrigation handle. Please do **not** take handle from its holder or splash off water!
- ☞ • Metal parts of the irrigation handle and of the double-barrelled hose will be hot! Please do not touch metal parts during the thermal cleaning procedure!
- During the thermal process for reducing bacterial count, the temperature display shows a 'd' on the first segment and the remaining time (in minutes) on the second segment.
- After cooling down to 37°C the unit will automatically switch into irrigation mode.

4.5.2 Irrigation mode

Purpose:

Rinsing of the auditory canal.

Properties:

- Temperature: 37°C, permanently set.
- Position of the adjusting lever: **High** (see fig. 6, page 26).

Activation:

- Automatically after completion of thermal cleaning procedure and taking out the handle or by key (Ⓞ, fig. 1).
- When the irrigation mode is not used (handle remains in holder) the Variotherm plus changes into the energy saving mode after 5 min.
- Reactivation of the irrigation mode by taking the handle out of the holder or activating a key on the temperature control unit.
- Repeated activation of the 37°C key effects switching-off of the heating system.
 - ↳ Water with feed temperature is available.

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Stimulation of the vestibule.

Properties:

- Temperature: according to pre-set cold or warm water stimulation level
- Position of the adjusting lever: **Low** (see fig. 6, page 26)
- Duration: as pre-set by means of timer.

Activation:

- Initial selection of the type of stimulation is done by actuating either warm-stimulation or cold-stimulation key (7, 8, fig. 1) (see section 4.1 for pre-setting of temperature).
 - ↳ Water with pre-set temperature value flows to handle.
- Wait for 20 seconds.
- Operate the "timer-start-key".
- Preparation for stimulation:
 - Operate unit 15 seconds in holding position, so that the jet can be positioned correctly in the auditory canal.

↳ During this period, the water is directed, via the handle, into the drain in order to guarantee an optimum temperature stability; with an optical indication by the flashing figure in temperature display.
Do not start stimulation by operating the release lever at the handle until you hear the audible signal!

- Carrying out the thermo-stimulation by actuating release lever at handle (9, fig. 4).
- At the end of the stimulation period a control signal for a recording unit is issued at the nystagmograph output.
- After completion of the stimulation the unit remains in a holding position for 15 seconds so that the handle may be removed from the auditory canal.
 - ↳ During this time, water flows into the drain; again an optical indication is made by flashing temperature value.
- Repeated actuation of the "timer-start-key" leads to stopping the timer operation.
- Second actuation of the currently active key effects deactivation of the pertaining stage.
 - ↳ Heating is switched off completely.
 - ↳ Stimulation with cold water (feed temperature).

4.5.4 Energy saving mode

Purpose:

- Reduction of energy consumption.

Activation:

- Actuation of key "heating on/off" (2, fig. 1).
 - ↳ Heating is switched off.
- After five minutes of non-use of the unit.

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5.0 Cleaning

5.1 General information on cleaning and disinfection

The jet connection which comes into contact with the patient must be disinfected after each use. The hose tips must be exchanged after each use.

The surfaces of the Variotherm plus resist most of the common surface disinfectants.

However, do not use

- disinfectants containing concentrated organic or anorganic acids or bases, since these may cause corrosion damage;
- disinfectants containing chloramides, phenol derivatives or anionic tensides, since these may cause stress cracks in the plastic material used for the housing of the unit.

For cleaning and disinfection, you may also use disinfectant sprays or disinfectant tissues.

- ☞ Switch off the unit by the main power switch before starting cleaning and disinfection.

Wipe the unit surface with a cloth moistened with a cleaning or disinfecting solution. Make sure that no liquid penetrates the unit. You may use all cleansers and disinfectants listed in section 5.3.

- ☞ Spilled liquid must immediately be dry wiped.
- ☞ Always observe the instructions for use by the manufacturer of the disinfectants, above all, concentration specifications.
- ☞ The above measures for cleaning and disinfecting do not replace the respective precautions for operating the unit.

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5.2 Recommended disinfectants

Recommended instrument disinfectants	Ingredients	(in 100 g)	Manufacturer
Korsolex basic (Application concentrate)	Glutaral (ethylene-dioxy)dimethanol tensides, salts, corrosion inhibitors	15.2 g 19.7 g	Bode Chemie, Hamburg
Sekusept aktiv (Application concentrate)	sodiumpercarbonate, phosphonates non-ionic tensides		Ecolab, Düsseldorf
Gigasept FF new (Application concentrate)	succinic acid dialdehyde dimethoxy tetrahydrofurane corrosion inhibitors	11.0 g 3.0 g	Schülke & Mayr, Norderstedt
Instruments, automatic			
Neodisher MediClean forte (Application concentrate)	NTA non-ionic tensides enzymes preserving agents	5-15 g < 5g	Dr. Weigert, Hamburg
Machine disinfection	Ingredients	(in 100 g)	Manufacturer
Dismozon pur (Application concentrate)	magnesium monoperoxyphthalate hexahydrate	80 g	Bode Chemie, Hamburg
Green & Clean SK (Application concentrate)	alkyl-dimethyl-benzyl-ammonium chloride dialkyl-dimethyl-ammonium chloride Alkyldimethylethylbenzylammoniumchlorid	<1 g	Metasys, Rum (Austria)
Perform	Pentakalium bis(peroxymonosulfate)- bis(sulfate)	45.0 g	Schülke & Mayr, Norderstedt

If using aldehyde-containing or amine-containing disinfectants at the same object, this may result in discolourations.

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
⁶⁾ Material bei 134°C formstabil

5.3 Cleaning method for handle with jet connection for water irrigation tip

The ATMOS appliances Variotherm and ATMOS Hygrotherm are used, according to the rules, with disposable hose tips (disposables). These tips must be disposed off after each application to the patient.

When using the hose tips ATMOS recommends a reprocessing in line with the below cleaning and disinfection plan.

	What	How			Recommendations	When				Who Qualified and trained staff who are familiar with reprocessing. (Please fill in the responsible person -> use a water-based overhead marker)
	Reusable parts	C Cleaning	D Disinfection	S Sterilisation		NAfter each procedure	Daily	Weekly	Monthly	
Ear irrigation / Thermal nystagmus stimulation										
	Ear irrigation bowl	X	X ^{2,4,5}		Manual or automatic cleaning and disinfection	X				
	Handle	X	X ³		Wipe cleaning and wipe disinfection		X			
	Jet connection	X	X ^{2,4,5,6}		Manual or automatic cleaning and disinfection		X			
	Splash protection	X	X ^{2,4,5}		Cleaning and disinfection (automatic or manual)		X			
	Hose tip (disposable)				Exchange after each application	X				
	Rinsing attachment	X	X ^{2,4,5}		Manual or automatic cleaning and disinfection	X				
	Hygiene filter				See operating instructions for hygiene filter				X	
	Rinsing lid with rinsing hose	X	X ^{2,4,5}		Manual or automatic cleaning and disinfection		X			
	Rinsing bottle	X	X ^{2,4,5,6}		Cleaning and disinfection; cleaning in the dishwasher with the glass care programme		X			

Recommended disinfectants

³⁾ Surface disinfection for coated surfaces:

- Green & Clean SK (ATMOS)
- Dismozon® pur (Bode Chemie)
- Kohrsolin® FF (Bode Chemie)
- Perform® (Schülke & Mayr)
- Terralin® Protect (Schülke & Mayr)

Other surfaces:

- Dismozon® pur (Bode Chemie)
- Kohrsolin® FF (Bode Chemie)
- Bacillocid® rasant (Bode Chemie)
- Mikrobac® forte (Bode Chemie)
- Perform® (Schülke & Mayr)
- Terralin® Protect (Schülke & Mayr)
- Surface disinfection FD 312 (Dürr Dental)
- Quick disinfection B 30 (Orochemie)

⁴⁾ Instruments - manual disinfection:

- Korsolex® AF (Bode Chemie)
- Korsolex® basic (Bode Chemie)
- Korsolex® plus (Bode Chemie)
- Korsolex® extra (Bode Chemie)
- neodisher® Septo MED (Dr. Weigert)
- neodisher® Septo 3000 (Dr. Weigert)
- Sekusept® PLUS (Ecolab)
- Sekusept® aktiv (Ecolab)
- Gigasept® Instru AF (Schülke & Mayr)
- Gigazyme® (Schülke & Mayr)
- Gigasept FF neu (Schülke & Mayr)

⁵⁾ Instruments - automatic disinfection:

- Dismoclean® 21 clean (Bode Chemie)
- Dismoclean® 24 Vario (Bode Chemie)
- Dismoclean® 28 alka one (Bode Chemie)
- Dismoclean® twin basic/twin zyme (Bode Chemie)
- neodisher® FA (Dr. Weigert)
- neodisher® FA forte (Dr. Weigert)
- neodisher® MediClean forte (Dr. Weigert)
- Thermosept® alka clean forte (Schülke & Mayr)
- Thermosept® RKN-zym (Schülke & Mayr)

Please see the manufacturer's instructions for concentration, contact time, temperature and the compatibility of materials.

Important information

Wipe cleaning and wipe disinfection: All surfaces have to be wiped with a clean (disposable) wipe which is damped with disinfectant solution; the entire surface has to be wiped thoroughly and may not be dried afterwards.

¹⁾ Please observe the manufacturer's operating instructions.

²⁾ Alternative to manual cleaning und disinfection: Wash-Disinfectant 78°C / 172°F

⁶⁾ Material dimensionally stable at 134°C

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6.0 Maintenance and service

Decalcification and changing of filter (see sections 6.4, 6.5) should be carried out immediately by the user when required.

Safety-technical inspection

6.1 Functional test

Test

- Switch unit several times on and off.
 - ↳ Correct function of mains pilot lamp and solenoid valve (clacking).
 - ↳ Automatic activation of the thermal process for reducing bacterial count.

- Shut-off main water tap and observe temperature indication.
 - ↳ After a few seconds error signal will appear "F1" (= no water).

- Testing the displays.
 - ↳ All optical and audible indicators issue their signal during the first second after switching on (in particular segmental test "8 8").

- Set all temperature setting switches on 40°C and test consecutively by touching the corresponding keys whether in all ranges 40°C is reached in ACTUAL temperature indication.
 - Spray off and check water temperature by hand.
 - Carry out the same test with temperatures 37°C and 28°C.
 - ↳ The particular combination of temperatures 40°C, 37°C and 28°C, as well as the test in all ranges ensures a comprehensive test of all contacts of the pre-setting switches, in ON- as well as in OFF-position, and also testing of all components involved in temperature control of the DESIRED temperature control system. Therefore, carrying out a test in one range only, or at only one temperature is not permissible.

- Set temperature at 49°C.
 - ↳ The electronic excess-temperature protection must respond (in the temperature display error code "F7" appears).

- Measure temperature of water splashed off at standard settings: 30°, 37° and 44°C.
 - ↳ The value may not deviate more than $\pm 0,5^\circ \pm 1/2$ digit, i.e. 1°C in total. Typically, however, values are better than $\pm 0,5^\circ\text{C}$.

This measurement is critical !

Measuring probes or thermometers which are too large may cause faulty measurements because of substantial heat removal, inertia, and because their temperature sensitivity is not localised on a point but spread over a larger surface. Use measuring instruments with an immersion probe, having a minimum precision of 0.3°C and a resolution of 0.1°C.

After setting the temperature, wait at least 1 minute before measuring to allow temperature stabilization of the double-barrelled hose and handle.

See also section 4.5.3!

6.2 Visual inspection

Test

- Check unit for external damage (damage to cables, hoses, filter and irrigation handle).
- Check hose system for external indications of leaks (drops, water spots, corrosion deposits etc.).
- Inspect unit, disconnected from the main power supply and with open housing, for internal damage, leaks, corrosion, and deposits (mineral scale). Close housing again and tighten screws. Re-connect main power cable.
- Check fuses for correct protection value and melting characteristics.
- Check sign plates and inscriptions for completeness and readability.

6.3 Electrical test

Test

- Measure **protective conductor resistance** according to VDE 0751 (max. 0.3 Ohm).
- Measure and compare with "first value measured" of **unit leakage current** according to VDE 0751 (measured during operation).
 - ↳ The "value measured first" may not be exceeded by more than 50%.
- Document and record the evaluation and tests according to VDE 0751.

6.4 Decalcification

- If the tap water in the area, where the unit is used, is known for mineral precipitation, preventive decalcification should be carried out:
 - Close water tap.
 - If applicable, dismantle filter when mounted vertically at water tap.
 - Unscrew screwed glass of filter.
 - Fill screwed glass with approx. 100 ml of decalcifier on acetic / citric, or lactic acid base (for inst. Citrosteril, by Fresenius).
 - Screw glass on again.
 - Install filter vertically again, if applicable.
 - Open water tap again.
 - Switch on the unit.
 - ↳ Decalcification takes place automatically in the thermal process for reducing bacterial count. After the unit has been in this mode for 5 minutes the water passages have been decalcified and cleaned, and also sufficiently flushed. The Variotherm plus is again fully ready for operation. When water passages are heavily calcified, a modified procedure may become necessary, with a prolonged reaction period for the decalcifier.
- Recommendation: Calcification safety device 502.0995.0.


Fig. 7. Water filter

- ❶ Filter body
- ❷ Filter cartridge
- ❸ Sealing ring
- ❹ Plexiglass cup

6.5 Replacement of filter

- When the filter cartridge is very dirty (dark colour) the filter cartridge is to be exchanged as follows:
 - Close water tap of supply.
 - Switch on unit to reduce pressure in the filter body.
 - Switch off unit after appearance of "F 1".
 - Unscrew Plexiglass cup (❹, fig. 7) from filter body (❶, fig. 7) by turning anti-clockwise.
 - Extract filter cartridge (❷, fig. 7) and replace by new one (Art.No. 502.0891.0).
- The filter element must not come into contact with contaminated objects in order to prevent the ingress of germs into the water system.

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

7.0 Trouble-shooting

 If, errors cannot be corrected inform the ATMOS service staff. Do not attempt to repair the unit yourself!

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8.0 Accessories and spare parts

8.1 Accessories

<u>Description</u>	REF
Calcification safety device.....	502.0995.0
Jet connection extra long, straight (110 mm).....	508.0429.0
Rinsing connection (Teflon) for stimulating the Vestibular organ	501.0355.0
Extension of supply hose G3/4a-G3/4i, L =1,5 m	501.0315.6
Waste water hose, special length, G1/4a-G1/4a (please indicate desired length)	502.0881.1
Drainage fittings as to DIN 1988	502.0880.0
Connecting cable for a nystagmograph	502.0850.0
Connecting cable for the equipotential compensation, L = 5m	008.0596.0

8.2 Spare parts

<u>Description</u>	REF
Variotherm plus , basic unit, 220-240 V AC, 50/60 Hz	502.0900.0
Variotherm plus , basic unit, 110-127 V AC, 50/60 Hz	502.0900.1
Irrigation handle 4+	502.0963.0
Double-barrelled hose.....	510.0412.0
Jet connection short, straight (80mm).....	508.0427.0
Splash protection	501.0331.0
Rubber bush for splash protection.....	501.0331.1
Throttling nozzle 0.5 mm (150 ml/min).....	502.0946.0
Throttling nozzle 0.7 mm (250 ml/min).....	502.0946.1
Three-holes gasket	501.0320.0
Hose tips for jet connection (30 pcs.).....	502.0844.0
G3/4a-G3/4i, L = 3 m	502.0768.1
Cap nut G3/4 (to be included in order).....	052.0035.1
Water filter complete	502.0890.0
Filter element for water filter	502.0891.0
Gasket (for G1/4a) 13x18x1.....	055.0018.0
G1/4a-G3/4a reducer (filter G1/4 to water hose G3/4).....	502.0724.1
Supply hose with connection G3/4i, L = 3 m.....	501.0315.7
Waste water hose with connection G1/4, L = 3 m.....	502.0882.0
Gasket (für G1/4i) 10x15x1	055.0003.0
Power cable with protective contact socket (Germany), L = 3m.....	507.0859.0
Power cable with British plug BS1363A (UK), L = 2,5m.....	008.0762.0
Fuse F 8 A (quick-acting) 250 V breaking capacity H (1500 A).....	008.0767.0
Fuse for special voltage (110-127VAC): F 15A (quick-acting) breaking capacity 10kA	008.0766.0
Front foil Variotherm plus	060.0357.0
Operating instructions.....	502.0900.i

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9.0 Specifications

Voltage range	220 - 240 V~ ± 10 % 50/60 Hz Special voltage: 100 - 127 V~ ± 10 % 50/60 Hz
Current input	6.8 A (220 - 240 V~); 14.8 A (100 - 127 V~)
Power cable	3 m
Power consumption	1500 W
Operating modes	Thermal process for reducing bacterial count; Irrigation mode; Stimulation mode; Heating off (energy saving)
Stimulation time	adjustable by timer from 1 up to 99 sec.
Timer indication	Accuracy ± 0.5 s ± ½ digit
Water temperature	20°C - 47°C
Temperature indication	Accuracy ± 0.6 °C ± ½ digit (can only be guarantee with unchanged ambient conditions)
Feed temperature of the water	+7...+19.5°C
Water quality	Drinking-water
Water flow in the irrigation mode (high flow)	450 ml/min: standard unit and variant 250/450 ml/min 400 ml/min: variant 400/400 ml/min 500 ml/min: variant 500/500 ml/min
Water flow in the stimulation mode (low flow)	150 ml/min: standard unit 150/450 ml/min (adjustable) 250 ml/min: variant 250/450 ml/min (adjustable) 400 ml/min: variant 400/400 ml/min (adjustable) 500 ml/min: variant 500/500 ml/min (adjustable)
Operating pressure	at least: +2000 hPa (2 bar); max: +6000 hPa (6 bar)
Connections	Power connection ; control output for a nystagmograph; equipotential equalization; water supply G 3/4" external thread; water drainage G 1/4" external thread; connection for double-barrelled hose
Operating time	Continuous operation
Fuses	quick-acting 8 A (220-240 V~); quick-acting 15 A (110-127 V~)
Protective earth conductor resistance	0,1 Ω max.
Earth leakage current	0.5 mA max.
Enclosure leakage current	0.1 mA max.
Patient leakage current	0.1 mA max.
Ambient conditions	
Transport/storage	-20...+50°C; 5...90 % humidity without condensing, air pressure 700...1060 hPa
Operation	+10...+35°C; 20...80 % humidity without condensing, air pressure 700...1060 hPa
Dimensions HxWxD	14.5 x 37 x 31.5 cm
Weight	approx. 5.5 kg
Protection class (EN 60601-1)	I
Degree of protection	Type B 
Protection category	IPX0
Classification acc. to Annex IX EEC directions 93/42/EEC	Ila
CE marking	CE 0124
Rules applied	EN 60601-1: 1990 + A1:1993 + A2:1995 EN 60601-1-2: 1993 (EMV / EMC)
UMDNS-Code	10-548
GMDN-Code	34891
REF	502.0900.0
Canadian Classification	
Device group	Ear, nose, throat
PNC	77ETP
Risk Class	2
Description	STIMULATOR, CALORIC (WATER)

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10.0 Disposal

- Packaging material, cardboard and/or PE foam, can be fully recycled or returned to your supplier.
- The Variotherm plus does not contain any hazardous materials.
- The material of the housing can be completely recycled.
- The component parts of the Variotherm plus must be correctly discarded and the materials should be carefully separated for recycling.
- Electronics circuit boards must be discarded using the appropriate recycling process.
- Used hose tips, which no longer can be disinfected, must be discarded into domestic waste immediately.

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This document is subject to change without notice.

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- i** ■ Medical electrical equipment is subject to special precautions with regard to EMC and must be installed according to following EMC notes.
- Portable and mobile HF communication facilities can influence medical electrical equipment.
- The use of other accessories, other converters and cables than stated may lead to an increased emission or a reduced interference immunity of the equipment or system.

11.1 Guidelines and Manufacturer's Declaration - Emissions

The ATMOS Variotherm plus is intended for use in the electromagnetic environment specified below. The customer or user of the ATMOS Variotherm plus should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The ATMOS Variotherm plus uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The ATMOS Variotherm plus is suitable for use in all establishments, including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	Class B	
Flicker IEC 61000-3-3	is conform	

- i** The device may not be used directly next to other devices or piled up with other devices. If operation next to or piled with other devices is necessary, please watch the device to check its intended operation in this arrangement.

11.2 Guidelines and Manufacturer's Declaration - Immunity

The ATMOS Variotherm plus is intended for use in the electromagnetic environment specified below. The customer or user of the ATMOS Variotherm plus should ensure that it is used in such an environment.

Immunity Test	IEC 60601-Test Level	Compliance Level	Electromagnetic Environment - Guidance
ESD IEC 61000-4-2	± 6 kV Contact ± 8 kV Air	± 6 kV Contact ± 8 kV Air	Floors should be wood, concrete, or ceramic tile. If floors are synthetic, the relative humidity should be at least 30%.
EFT IEC 61000-4-4	± 2 kV Mains ± 1 kV I/Os	± 2 kV Mains Inapplicable	Mains power quality should be that of a typical commercial or hospital environment.
Surges IEC 61000-4-5	± 1 kV Differential ± 1 kV Common	± 2 kV Differential ± 1 kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Power Frequency 50/60 Hz Magnetic field IEC 61000-4-8	3 A/m	Inapplicable	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

Immunity Test	IEC 60601-Test Level	Compliance Level	Electromagnetic Environment - Guidance
Voltage Dips / Dropout IEC 61000-4-11	<p>< 5 % U_T (> 95 % Dip of the U_T) for 0.5 Cycle</p> <p>40 % U_T (60% Dip of the U_T) for 5 Cycles</p> <p>70% U_T (30 % Dip of the U_T) for 25 Cycles</p> <p>< 5 % U_T (>95 % Dip of the U_T) for 5 s</p>	<p>< 5 % U_T (> 95 % Einbruch der U_T) für 0,5 Periode</p> <p>40 % U_T (60% Einbruch der U_T) für 5 Perioden</p> <p>70% U_T (30 % Einbruch der U_T) für 25 Perioden</p> <p>< 5 % U_T (>95 % Einbruch der U_T) für 5 s</p>	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ATMOS Variotherm plus demands continued function even in case of interruptions of the energy supply, it is recommended to supply the ATMOS Variotherm plus from an uninterruptible current supply or a battery.
NOTE U _T is the mains alternating current prior to application of the test levels.			

11.3 Guidelines and Manufacturer's Declaration - Immunity

The ATMOS Variotherm plus is intended for use in the electromagnetic environment specified below. The customer or user of the ATMOS Variotherm plus should ensure that it is used in such an environment.

Immunity Test	IEC 60601-Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	3 V	Portable and mobile communications equipment should be separated from the ATMOS Variotherm plus incl. the cables by no less than the distances calculated/listed below. Recommended distances: $d = (3,5 / \sqrt{P}) * \sqrt{(P)}$ $d = (3,5 / \sqrt{E1}) * \sqrt{(P)}$ 80-800 MHz $d = (7 / \sqrt{E1}) * \sqrt{(P)}$ 0,8-2,5 GHz where „P“ is the max. power in watts (W) and D is the recommended separation distance in meters (m). Field strengths from fixed transmitters, as determined by an electromagnetic site (a) survey, should be less than the compliance level (b). Interference may occur in the vicinity of equipment containing following symbol.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1	With 80 MHz and 800 MHz the higher frequency range applies.
NOTE 2	These guidelines might not be applicable in any case. The propagation of electromagnetic sizes is influenced by absorptions and reflections of buildings, objects and people.
a	The field strength of stationary transmitters, such as base stations of cellular phones and mobile terrain radio equipment, amateur radio transmitters, cbm broadcast and TV stations cannot be predestined exactly. To determine the electromagnetic environment in regard to stationary transmitters, a study of the location is to be considered. If the measured field strength at the location where the ATMOS Variotherm plus is used exceeds the above compliance level, the ATMOS Variotherm plus is to be observed to verify the intended use. If abnormal performance characteristics are noted, additional measures might be necessary, e. g. a changed arrangement or another location for the device.
b	Within the frequency range of 150 kHz to 80 MHz the field strength is to be below 3 V/m.

11.4 Recommended separations between portable and mobile RF Communications equipment and the ATMOS Variotherm plus

The ATMOS Variotherm plus is intended for use in electromagnetic environment in which radiated disturbances are controlled. The customer or user of the ATMOS Variotherm plus can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications equipment and the ATMOS Variotherm plus as stated below, according to the maximum output power of the communications equipment.

Nominal output of the transmitter W	Separation distance, depending on transmit-frequency m		
	150 kHz bis 80 MHz $d = [3,5 / 3] \sqrt{P}$	80 MHz bis 800 MHz $d = [3,5 / 3] \sqrt{P}$	800 MHz bis 2,5 GHz $d = [7,0 / 3] \sqrt{P}$
0.01	0.12	0.12	0.24
0.1	0.37	0.37	0.74
1	1.2	1.2	2.4
10	3.69	3.69	7.38
100	11.66	11.66	23.32

For transmitters for which the maximum nominal output isn't indicated in the above table, the recommended separation distance d in meters (m) can be determined using the equation belonging to the respective column whereas P is the maximum nominal output of the transmitter in watts (W) acc. to manufacturer's specification.

NOTE 1 With 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2 These guidelines are not applicable in any case. The propagation of electromagnetic sizes is influenced by absorptions and reflections of buildings, objects and people.

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1. General:

Our General Standard Terms and Conditions apply exclusively. Client's terms and conditions which are contrary to or deviate from our General Standard Terms and Conditions are not recognised unless their validity is explicitly confirmed in writing. Our General Standard Terms and Conditions also apply even if we deliver to clients without reservation, in the knowledge of the client's contrary terms and conditions. Our General Standard Terms and Conditions also apply to all future business with that client.

2. Proposal - Order Confirmation

Our proposals are subject to change without notice unless otherwise stated in our order confirmation. Each order is only accepted by us following our written order confirmation.

3. Orders

Every order requires an exact description of all of our product's details. We assume no liability for errors and damage caused by inaccurate or incomplete ordering details.

4. Prices

Unless otherwise stated in the order confirmation, our prices in the order confirmation are ex factory prices and exclude packaging and value added tax. Packaging is charged separately at cost price in the invoice. Value added tax is charged separately in the invoice according to the legal rate on the invoice date. We reserve the right to change prices appropriately should price reductions or increases, especially due to wage settlements, changes in the price of materials or currency fluctuations, be incurred. Proof of such changes will be provided for the client on request.

5. Payment Conditions - Balancing

Unless otherwise stated in the order confirmation, our invoices are payable with a 3% discount within 10 days (except for repair and assembly services) or within 21 days from the invoice date net cash; money receipts is decisive for complying with this term. We are entitled to charge interest after the due date at a rate 2% above the relevant basic interest rate of the German Federal Bank. Should the client have payment arrears, we are entitled to charge interest on arrears at a rate 5% above the relevant basic interest rate of the German Federal Bank. Should we be able to prove higher damages due to arrears, we are also entitled to claim these. The client only has the right to balance invoices against its own claims should such claims be confirmed in a court of law or recognised by us. The client does not have the right of retention due to disputed counterclaims.

6. Delivery Periods

Fulfillment of our delivery duties requires the punctual and proper fulfillment of the client's duties. The right to defense on the grounds of an unfulfilled contract is reserved. Should the client default in accepting the goods delivery or breach other cooperation duties, we are entitled either to withdraw from the contract or claim compensation for any increased costs incurred up to that time without setting a further deadline. The right to make further claims is reserved. Furthermore, in such cases, the risk of coincidental destruction or a coincidental deterioration in the quality of the delivered goods is transferred to the client in the case of default in accepting such goods or payment arrears. Acts of God or stoppages (due to insufficient supplies of material, industrial disputes etc.) entitle us either to demand an appropriate extension of delivery periods or to partly or entirely dissolve the delivery contract. This does not give the client the right to claim damages. We have fulfilled delivery periods if the delivery goods have left our factory or the client has been informed of the goods' readiness for delivery within such delivery periods. Delivery periods stipulated by the client are not recognised by us unless they form part of our order confirmation. We adhere to legal terms and conditions in cases where, as a result of an undue delay in the delivery for which we are liable, the client is entitled to claim that his interests in a continued fulfillment of the contract have ceased. We also adhere to legal terms and conditions should a delay in delivery be caused by deliberate or grossly negligent action by us or our representatives for which we are responsible. We are also responsible for such actions by our representatives or agents. Should the delivery delay not be caused by our deliberate infringement of contractual duties for which we are responsible, our liability is limited to damage which is regarded as typical for that case. We are liable according to the legal terms and conditions if and in so far as the delivery delay for which we are responsible is caused by an infringement of a substantial contractual duty. In such cases, our liability is also limited to damage which is regarded as typical for that

case. Should the delivery delay be caused by a culpable infringement of non-substantial contractual duties, our client is also entitled to claim a one-off damage compensation worth 3 percentage points of the delivery value of the goods for each week's delay, up to a maximum which is no higher than 15 percentage points of the delivery value of the goods

7. Delivery - Familiarisation

In the case of the delivery of devices for the medico-technical industry which require assembly and/or familiarisation for the final customer using specialist trade personnel (such as Ear, Nose and Throat Apparatus and Suction Units), we reserve the right to deliver the goods exclusively to the relevant specialist traders. Should the trader not carry out assembly and/or familiarisation for the final customer, this is carried out by us. In such cases, we reserve the right to charge the client for the additionally created costs. Our specialist traders operate a recording system so that, if necessary, our products can be traced to the final customer. The specialist trader undertakes to immediately report to us all events and risks which must be reported in connection with our products.

8. Passage of Risk - Packaging

Unless otherwise stated in our order confirmation, delivery is agreed ex factory. The risk of the goods' damage or loss is therefore transferred to the client as soon as the goods leave the factory or the client is in default of acceptance of the goods. This also applies to cases where we confirm prepaid carriage. Transport packaging and all other packaging according to the packaging regulations is not returnable. Our client is responsible for disposing the packaging at its own cost. Our deliveries are insured by us at the client's expense unless explicitly otherwise agreed. No insurance is arranged in the case of goods which are collected by our clients. In the case of transport damage, claims are only handled if the client receives confirmation of any damage, reduced weight or loss by the shipping company before accepting the delivery.

9. Warranty

The client is responsible for examining the delivered goods immediately after receiving them to determine any eventual deficiencies or delivery errors, and to report these immediately. Should the client fulfil this examining and reporting responsibility, and should payment conditions be fulfilled, we shall be liable to the client within the scope of legal regulations. Our period of warranty shall in all cases be two years. Our client can make use of the warranty as follows, so long as he can provide first buyer proof (in the form of an invoice or delivery note) and provided that the product still has the original, unchanged serial number:

- a. We choose whether to fulfil our guarantee by providing repair services free of charge - either on the client's premises or in our factory - or replacing the product. We can also provide these guarantee services through an authorised company;
- b. Should a product be returned to us, the client agrees to send the product in its original or similar packaging, offering the same protection as the original packaging, to our address or any address notified by us.
- c. Our guarantee ceases to apply if changes of any kind have been made to our product, unless such changes have been made by us or a company authorised by us, or have been previously agreed upon in writing by us. Our guarantee also ceases to apply if third parties have carried out repairs to our products or replaced parts thereof. This applies regardless of the fact whether these measures individually or collectively led to a deficiency of the product;
- d. We accept no responsibility for damage defects caused by
 - operational wear and tear;
 - incorrect installation or incorrect or insufficient maintenance;
 - incorrect operation of the product (in contradiction to the operating instructions delivered with the product);
 - improper use or operating faults;
 - inappropriate or negligent handling and care, especially with respect to dirt, lime, suction of fluids, inappropriate cleaning and sterilisation;
 - using accessories and/or replacement parts which are not explicitly approved;
 - incorrect assembly and/or initial operation by the client or third parties;
 - the client's negligence in handling the product;
 - unacceptable operating conditions, such as humidity, temperatures, the power supply, vibrations.
 - accidents, acts of God, especially lightning, water, fire, public unrest and insufficient ventilation. We are not liable for damage to other objects apart from our product itself, except in the case of any deliberate or grossly negligent actions by us or our representatives or agents. Should no deliberate breach of contract be claimed, our liability

is limited to damage which is regarded as typical for that case. This also applies in the case of our culpable infringement of substantial contractual duties. The indispensable conditions of German Liability Law remain unaffected thereby.

- For second-hand equipment, the period of warranty shall be reduced to a period of twelve months.

10. Reservation of Ownership

We retain ownership of our goods until the receipt of all payments arising from the business relationship, including all demands arising from installation orders, subsequent orders, repairs, accessory deliveries and replacement orders. Should we have agreed upon payment on the basis of cheque and bill transactions, the ownership reservation applies until the cheque received by us has been paid in, and does not expire through our credit upon receiving the client's cheque. In the case of a breach of contract by the client, especially payment arrears, we are entitled to repossess our goods. Repossession of our goods represents a withdrawal from the contract, unless explicitly declared in writing by us. We have the right to utilise the product after its repossession, whilst the income from such use is balanced against the client's arrears, after deducting appropriate utilisation costs. The client is responsible for handling the goods with care. Should maintenance and inspection work be necessary, the client must carry these out punctually at his own cost. Our client is entitled to sell the goods he has bought from us in a proper sale transaction. However, he must immediately assign all outstanding claims to the value of the final invoice sum (including value added tax) of our claims to his customers or third parties. The client is entitled to collect this claim even after such assignment. Our right to collect the claim ourselves remains unaffected thereby. We undertake to release the securities to which we are entitled if requested to do so by the client should the realisable value of the our securities be more than 10 percentage points higher than the outstanding claims. We reserve the right to choose the securities to be released.

11. Plans and Illustrations

We retain ownership of and copyrights to all plans, illustrations, calculations and other documents which are attached to our proposals. The client must receive explicit written permission before passing these on to third parties. Imitating our legally patented products is forbidden and will be prosecuted.

12. Jurisdiction and Place of Performance

Our central office is the place of performance for all disputes in connection with these General Standard Terms and Conditions and the contracts closed with clients under them. This jurisdiction excludes other jurisdiction relating to persons or subject-matter. Furthermore, our client is not entitled to bring charges against us in another court should he file counter-charges, carry out counterbalancing or declare retention. We, however, are entitled to bring charges against our client at their general place of jurisdiction or at another relevant court recognised by German or foreign law. Unless otherwise stated in the order confirmation, our central office is the place of performance.

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