MedizinTechnik
... for a better life


Deutsch

English

Français

Gebrauchsanweisung
Operating Instructions
Mode d'emploi

Varioair 3

Inhalt

Seite

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Table of contents

	Page
1.0 Introduction	
1.1 Notes on operating instructions.....	5
1.2 Intended use	7
1.3 Function	7
1.4 Explanation of symbols	9
2.0 Safety advice	11
3.0 Installation and start-up	
3.1 Standard delivery	13
3.2 Illustrations	14 - 15
3.3 Connections	17
3.3.1 Electrical connection	17
3.3.2 Connecting a nystagmograph	17
3.3.3 Connecting a potential balance.....	17
3.3.4 Connection of the handle	17
3.3.5 Air inlet	17
3.4 Start-up	17
4.0 Operation	
4.1 Adjusting temperatures	19
4.2 Selecting temperature levels.....	19
4.3 Adjusting stimulation time	19
4.4 Description of operating modes	21
4.4.1 Preparation mode.....	21
4.4.2 Stimulation mode	23
4.4.3 Standby mode	23
5.0 Cleaning	
5.1 General information on cleaning and disinfection	25
5.2 Recommended instrument disinfectants	27
5.3 Recommended surface disinfectants	27
6.0 Maintenance and service	
6.1 Functional test.....	31
6.2 Visual inspection	32
6.3 Electrical test.....	32
6.4 Replacing of filter	32
7.0 Trouble shooting	36
8.0 Accessories and spare parts	
8.1 Accessories.....	39
8.2 Spare parts.....	39
9.0 Specifications	42
10.0 Disposal	44

Notes on EMC

General Standard Terms and Conditions

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1.0 Introduction

1.1 Notes on operating instructions

- The operating instructions contain important information on how to operate the Varioair 3 safely, correctly and effectively. The operating instructions are intended for new operating personnel and also should be continued to be used as a reference manual. Reading and understanding the operating instructions will assist the user in avoiding risks, reducing repair costs, lowering down-time, increasing reliability and extending the service-life of the equipment. **The operating instructions must always be kept available and near the equipment.** Prior to first use, please read the chapter "Safety advice", to be prepared for any possible dangerous situations.

The basic principles are:

Judicious and careful work provides the best protection against accidents!

Operational safety and readiness for use depend not only on your capabilities, but also on care and maintenance given to the Varioair 3. Regular cleaning and service work are a must. Major maintenance and repair work should be carried out only by expert personnel authorised by ATMOS. For repairs, only use ATMOS spare parts to preserve the equipment warranty and ensure operational safety, readiness for work and to maintain the value of the equipment.

- The product Varioair 3 bears CE Marking CE-0124 according to the EEC guideline of the council for medical products 93/42/EEC and meets the basic requirements of this guideline.
- The quality management system applied at ATMOS has been certified according to international standards EN ISO 9001 and EN ISO 13485.
- Reprints, and text extracts are allowed only with written permission from ATMOS.

Abbreviations / symbols in these operating instructions:

- Indicating a list
 - Subdivision of a list/activity

The recommended sequence must be followed in each case!

 Indicating particularly important advice!

 Describing the effect of an activity.

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1.2 Intended use

- Compact unit for stimulation of the vestibule. Distinguished by its easy operation, it can be used in hospital and in a practice.
- This unit may not be used if the auditory canal is inflamed or contaminated.
- This unit is designed for short-time use on patients.

1.3 Function

- After activating the main switch, the optical readout displays tested.
- The unit automatically switches into the standby mode in which the heating and the pump are switched off.
- The unit can be changed into the stimulation mode for stimulating the vestibule. The Varioair 3 is equipped with a timer for preselecting the stimulation time.

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1.4 Explanation of symbols

-  Pay attention to operating instructions !
as to ISO /7000/0434 DIN 30600/1008
IEC 348
-  Type B equipment as to IEC 417
-  Fuse
as to IEC 417/5016, DIN 30600/0186
- °C** Temperature in degree centigrade
- S** Timer adjustment in seconds
-  Start
-  Stop
-  Timer
-  Cold stimulation level
-  Warm stimulation level
-  Heating and air flow ON
-  Heating and air flow OFF (standby)
-  Control output for connecting a nystagmograph
(graphical recorder as to DIN 30600,
IEC 417 5192)
-  Equipotentiality
DIN 30600 495, ISO 417 5021
-  Air filter
DIN 24300

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2.0 Safety advice

- The Varioair 3 is produced according to IEC 601/EN 60601 and listed in the following classes:
 - VDE Class of protection 1
 - Class IIa (EEC 93/42).
- The equipment must only be connected to a power outlet with a correctly installed ground/earthing contact according to local and state regulations.
- Regulations for users in the European Economic Area: The Varioair 3 must be subjected annually to a safety-technical inspection according to the regulations in these operating instructions.
- The Varioair 3 may be used in supervised operation by qualified personnel that have been authorised by ATMOS and been trained for operating the equipment (IEC 601-1/ EN 60601-1).
- The mains voltage specified on the type plate must match the power supply system.
- Prior to every use of the equipment, verify that the device is technically safe and in proper working condition. **Damaged leads** must be replaced immediately!
- Correct configuration in assembly of country-specific connections

green/yellow:	protective conductor (PE)
blue:	neutral conductor (N)
black resp. brown:	phase (L)
- The control panel must be visible and in reach of the operator. Ensure the equipment is stable prior to use.
- 
 ● Prior to application, the air temperature must be checked by the user!
- Switch off the main power switch after using the equipment.
- The Varioair 3 should only be operated in rooms used for medical purposes, and never in areas subject to explosion hazards.
- All additional equipment, which is connected to the analog and digital interfaces of the unit, must meet the requirements of relevant EN specifications (for inst. EN 60950 for data processing equipment and EN 60601 for electrical medical appliances). In addition, configurations must satisfy system specification EN 60601-1-1. When additional equipment is connected to the signal input or signal output section on the unit, the person making the connection is deemed "a system configuration operator" and as such is responsible for meeting the requirements of system specification EN 60601-1-1. For answers to additional questions, please contact your local specialist supplier or ATMOS technical service.
- ATMOS is not liable for personal injury and damage to property if
 - original ATMOS parts are not used,
 - the operating instructions are not completely followed,
 - assembly, new settings, alterations, extensions and repair have been carried out by personnel not authorised by ATMOS.
- The validity for the enclosed certificate of conformity expires if the customer or a third party manipulates the unit, e.g. modifications of any kind, by installing non-authorized accessories, removing warning or information labels as well as using the unit for inappropriate applications.
- Please note:
A medical insulating transformer with earth leakage monitor or any similar safety system acc. to EN 60 601-1 is required, if several devices are connected over one common power supply. The transformer must correspond to the power consumption of all the devices to be connected.

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3.1 Standard delivery

Quantity

1	Varioair 3 Basic unit
1	Power cable
1	Handle with tube
1	Hose tips (30 pcs.)
1	Operating instructions

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D

Fig. 1. Front view

- ① Main switch
- ② Key switch f. heating/air flow ON/OFF (standby)
- ③ Coding switch for warm-stimulation level
- ④ Coding switch for cold-stimulation level
- ⑤ Coding switch for stimulation time
- ⑥ Key switch for selecting the warm-stimulation level (e.g. 44°C)
- ⑦ Key switch for selecting the cold-stimulation level (e.g. 30°C)
- ⑧ Key switch for start/stop of the stimulation
- ⑨ Temperature display (two-figure number, resolution 1°C), current value indication
- ⑩ Display of stimulation time (two-figure number, resolution 1s)
- ⑪ Support for handle
- ⑫ Connection for handle

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D

Fig. 2. Rear view

- ❶ Control output for a nystagmograph
- ❷ Unit plug with fuse compartment
- ❸ Equipotentiality connection
- ❹ Air inlet filter

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D

Fig. 3. Handle

- ❶ LED for indicating the stimulating process
- ❷ Timer Start/Stop key

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3.3 Connections

3.3.1 Electrical connection

- Connect the power cable to the socket on the unit (②, fig. 2).
- Insert the power plug in a correctly installed socket with earthing contact.

3.3.2 Connecting a nystagmograph

- When controlling an ENG (electro-nystagmograph) or a CNG (computer-nystagmograph) at output (①, fig. 2), please connect only recording equipment approved by ATMOS. Connecting cable available from ATMOS (see section 8.0).

3.3.3 Equipotentiality connection

- Connection for potential compensation (③, fig. 2). Connecting cable available from ATMOS (see section 8.0).

3.3.4 Connection of the handle

- Only the handle which is intended for this purpose must be used:
 - Press the special connector plug slightly onto the socket in the unit.
 - Fasten it to the housing of the unit by turning the holding screws to the right.

 Do not buckle the air tube !

3.3.5 Air inlet

- The ambient air is sucked into the unit through the filter.
- Screw connector for the air filter.

3.4 Start-up

- Insert the handle in its holder; the air outlet must point to reverse side of the unit.
- Switch on the unit (①, fig. 1).
- Automatic display test with digital numbers "88" and audible signal.
- Automatic change into the standby mode.

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- Mittels Codierschalter (⑤, Bild 1).

4.0 Operation

- Use the main switch to start the unit.
- The unit is ready for operation after 1 second.

4.1 Adjusting temperatures

- Two variable temperature levels (20°C - 47°C) (48°C + 49°C only for testing purposes).

☞ The lowest achievable stimulation temperature is approx. 2°C over the ambient temperature.

- Temperate setting by coding switches (3, 4, fig. 1).
 - left switch: for adjusting the "ten" partition
 - right switch: for adjusting the "one" partition
 - ☞ lower keys (+):
Temperature increase
 - ☞ upper keys (-):
Temperature decrease
- Standard settings:
 - Level for cold stimulation: 30°C
 - Level for warm stimulation: 44°C

4.2 Selecting temperature levels

- For selecting the desired temperature level, use the respective key (6, 7, fig. 1).
 - ☞ Display of the active level by illuminated diodes.
 - ☞ Display of the air temperature (current value) in °C.
- For switching off the heating system press the respective key (6, 7, fig. 1) of the active temperature level.
 - ☞ Light of LED of temperature level disappears.
 - ☞ Air temperature is shown (actual value) in °C.

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4.4.1 Preparation mode

Purpose:

The temperature set by the user is adjustable.

Properties:

- Temperature: corresponds to the preselected cold or warm stimulation level
- Air flow: 5.0 l/min.

Activation:

- By operating one of the key switches for selecting a temperature level (6, 7, fig. 1) or by operating the key switch for heating/air flow ON/OFF (2, fig. 1) in the standby mode
- If the active temperature push button is repeatedly pressed, the heating is switched off.
 - ↳ Air with a temperature that almost corresponds to the ambient air temperature is available.

Deactivation:

- Operating the key switch heating/air flow ON/OFF (2, fig. 1) switches the unit to the standby mode.
- Automatically switches to the standby mode when the unit is not used within a 3 minute period.

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Stimulation of the vestibule.

Properties:

- Temperature: Set according to pre-set cold or warm stimulation level
- Air flow: 5 l/min
- Duration: As pre-set by the timer.

Activation:

- Initially select the type of stimulation by actuating either warm-stimulation or cold-stimulation key (6, 7, fig. 1) (see section 4.1 for pre-setting of temperature).
- Operate the "timer-start-key" on the unit or on the handle.
- Preparation for stimulation:
 - As long as the push button is pressed, the pump remains switched off to allow the jet connection to be positioned in the auditory canal.
- When the push button is released, the thermal stimulation is performed for the time set by the user (5, fig. 1).
 - ↳ LED on handle (1, fig. 3) lights up during thermo-stimulation.
- At the end of the stimulation period a control signal for a recording unit is issued at the nystagmograph output.
- After completion of stimulation the pump is switched off.
- Repeated actuation of the "timer-start-key" during the stimulation leads to stopping timer operation.
- Second actuation of the currently active key effects deactivation of the pertaining stage.
 - ↳ Heating is switched off.
 - ↳ Stimulation with near-ambient air temperature.

4.4.3 Standby mode

Purpose:

- Reduction of energy consumption.
- Reduction of noise level.

Activation:

- Actuation of key "heating/pump on/off" (2, fig. 1).
 - ↳ Heating is switched off.
 - ↳ Pump is switched off after 2 s.
- Automatically switches to the standby mode after each stimulation process.
- Automatically switches to the standby mode when the unit is not used within a 3 minute period.

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5.0 Cleaning

5.1 General information on cleaning and disinfection

After use, all parts which come into contact with the patient (jet connection, hose tips) must be removed and disinfected resp. autoclaved! You may use all cleansers and disinfectants listed in section 5.2 can be used.

The handle is not autoclavable!

The surfaces of the Varioair 3 resist most of the common surface disinfectants.

However, do not use

- disinfectants containing concentrated organic or anorganic acids or bases, since these may cause corrosion damages;
- disinfectants containing chloramides, phenol derivatives or anionic tensides, since these may cause stress cracks in the plastic material used for the housing of the unit.

For cleaning and disinfection, you may also use disinfectant sprays or disinfectant tissues.

 Switch off the unit using the main power switch before starting to clean and disinfect.

Wipe the unit surface with a cloth moistened with a cleaning or disinfecting solution. Make sure that no liquid penetrates the unit. You may use all cleansers and disinfectants listed in section 5.3.

 Always observe the instructions for use by the manufacturer of the disinfectants, including all concentration specifications.

 The above measures for cleaning and sterilizing do not replace the precautions for operating the unit.

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6.0 Maintenance and service

Change of filter (see section 6.4) should be carried out immediately by the user when required.

Safety-technical inspection

6.1 Functional test

Test

- Switch unit several times on and off.
 - ↳ Correct function of mains pilot lamp.
 - ↳ Automatic activation of displays test.
- Testing the displays.
 - ↳ All optical and audible indicators issue their signal during the first second after switching on (in particular segmental test "8 8").
- Set all temperature setting switches on 40°C and test consecutively by touching the corresponding keys whether in all ranges 40°C is reached in ACTUAL temperature indication.
 - Check air temperature.
 - Carry out the same test with temperatures 37°C and 28°C.
 - ↳ The particular combination of temperatures 40°C, 37°C and 28°C, as well as the test in all ranges ensures a comprehensive test of all contacts of the pre-setting switches, in ON- as well as in OFF-position, and also testing of all components involved in temperature control of the DESIRED temperature control system. Therefore, carrying out a test in one range only, or at only one temperature is not permissible.
- Set temperature at 49°C and close the air exhaust opening on the handle.
 - ↳ The electronic excess-temperature protection must respond after approx. 100-120 sec. (in the temperature display error code "F7" appears).
- Measure temperature of air at standard settings: 30° and 44°C.
 - ↳ The value may not deviate more than $\pm 0,5^\circ \pm 1/2$ digit, i.e. 1°C in total. Typically, however, values are better than $\pm 0,5^\circ$.

This measurement is critical !

Measuring probes or thermometers which are too large may cause faulty measurements because of substantial heat removal, inertia, and because their temperature sensitivity is not localised on a point but spread over a larger surface. Use measuring instruments with an immersion probe, having a minimum precision of 0.3°C and a resolution of 0.1°C.

After setting the temperature, wait at least 1 minute before measuring to allow temperature stabilization of the double-barrelled hose and handle.

After setting the temperature 1 minute at least should be waited in order that the probe attains correct temperature.

6.2 Visual inspection

Test

- Check the unit for external damage (damage to cables, hoses, filter and handle).
- Check filter for dirt or a defect.
- Inspect unit, disconnected from the main power supply and with open housing, for internal damage, and leaks. Close housing again and tighten screws. Re-connect the main power cable.
- Check fuses for correct protection value and melting characteristics.
- Check sign plates and inscriptions for completeness and readability.

6.3 Electrical test

Test

- Measure protective conductor resistance according to VDE 0751 (max. 0.3 Ohm).
- Measure and compare with "first value measured" of unit leakage current according to VDE 0751 (measured during operation).
 - ↳ The "value measured first" may not be exceeded by more than 50%.
- Document and record the evaluation and tests according to VDE 0751.

6.4 Replacing the filter

- The air filter (④, fig. 2) must be replaced once per year as follows:
 - Remove the filter by unscrewing it from the holder (turn the filter to the left).
 - Mount the new filter by screwing it into the holder (turn the filter to the right).

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

7.0 Trouble-shooting

If errors cannot be corrected with the assistance of the trouble-shooting list, please inform the ATMOS service staff. Do not attempt to repair the unit yourself!

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8.0 Accessories and spare parts

8.1 Accessories

<u>Description</u>	<u>REF</u>
Air inlet filter	000.0545.0
Power cable	507.0859.0
Handle, complete.....	502.1035.0
Hose tips for jet (30 pcs.).....	502.0844.0
Operating instructions.....	502.1100.i

8.2 Spare parts

<u>Description</u>	<u>REF</u>
Unit basis, to be plugged in.....	000.0639.0
Fuse T800mA H / 250V 5x20mm.....	008.0081.0
Air inlet filter	000.0545.0
Power cable	507.0859.0
Handle, complete.....	502.1035.0
Jet connection.....	502.1045.0
Hose tips for jet (30 pcs.).....	502.0844.0

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9.0 Specifications

Voltage range	100 - 240 V~ ± 5 %; 50/60 Hz
Current input	0.75 A max.
Power cable	3 m
Power consumption	85 W max.
Operating modes	Preparation mode; Stimulation mode (at temperature preselected for the cold stimulation level resp. warm stimulation level); Heating off and no air flow (economy mode, standby mode)
Stimulation time	can be preselected by timer from 1 up to 99 sec.
Timer indication	Accuracy ± 0.5 s ± ½ digit
Air temperature	20°C - 47°C
Lowest temperature	approx. 2°C above room temperature
Temperature indication	Accuracy ± 0.5 °C ± ½ digit
Temperature deviation	< ± 1°C
Air flow	5.0 l/min ± 10 %
Connections	Power connection; Control output for a nystagmograph; Equipotential equalization; Connection for the handle; Air inlet
Operating time	Short-time operation: 1. Automatic shut-off after completion of stimulation. 2. Automatic shut-off after 3 minutes.
Fuses	T 0,8 A
Protective earth conductor resistance	0.1 Ω max.
Earth leakage current	0.5 mA max.
Enclosure leakage current	0.1 mA max.
Patient leakage current	0.1 mA max.
Ambient conditions	
Transport/storage	-20...+50°C; 5...90 % humidity without condensing, air pressure 700...1060 hPa
Operation	+10...+35°C; 20...80 % humidity without condensing, air pressure 700...1060 hPa
Dimensions HxWxD	14.5 x 37 x 32 cm
Weight	approx. 3.7 kg
Protection class (EN 60601-1)	I
Degree of protection	Type B 
Protection category	IPX0
Classification acc. to Annex IX EEC directions 93/42/EEC	Ila
CE marking	CE 0124
Rules applied	EN 60601-1: 1990 + A1:1993 + A2:1995 EN 60601-1-2: 1993 (EMV / EMC)
UMDNS-Code	10-548
REF	502.1100.0
Canadian Classification	
Device group	Ear, nose, throat
PNC	77ETP
Risk Class	2
Description	STIMULATOR, CALORIC (WATER)

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10.0 Disposal

- Packaging material, cardboard and/or PE foam, can be fully recycled or returned to your supplier.
- The Varioair 3 does not contain any hazardous materials.
- The material of the housing can be completely recycled.
- The component parts of the Varioair 3 must be correctly discarded and the materials should be carefully separated for recycling.
- Electronics circuit boards must be discarded using the appropriate recycling process.
- Used hose tips, which no longer can be disinfected, must be discarded into domestic waste immediately.

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11.0 Notes on EMC

UK

-  ■ Medical electrical equipment is subject to special precautions with regard to EMC and must be installed according to following EMC notes.
- Portable and mobile HF communication facilities can influence medical electrical equipment.
- The use of other accessories, other converters and cables than stated may lead to an increased emission or a reduced interference immunity of the equipment or system.

11.1 Guidelines and Manufacturer's Declaration - Emissions

The ENT-camera Varioair 3 is intended for use in the electromagnetic environment specified below. The customer or user of the Varioair 3 should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The ATMOS Varioair 3 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The ATMOS Varioair 3 is suitable for use in all establishments, including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	Class A	
Flicker IEC 61000-3-3	Match	

-  The device may not be used directly next to other devices or piled up with other devices. If operation next to or piled with other devices is necessary, please watch the device to check its intended operation in this arrangement.

11.2 Guidelines and Manufacturer's Declaration - Immunity

The ENT-camera Varioair 3 is intended for use in the electromagnetic environment specified below. The customer or user of the Varioair 3 should ensure that it is used in such an environment.

Immunity Test	IEC 60601-Test Level	Compliance Level	Electromagnetic Environment - Guidance
ESD IEC 61000-4-2		± 6 kV Contact ± 8 kV Air	Floors should be wood, concrete, or ceramic tile. If floors are synthetic, the relative humidity should be at least 30%.
EFT IEC 61000-4-4		± 2 kV Mains ± 1 kV I/Os	Mains power quality should be that of a typical commercial or hospital environment.
Surges IEC 61000-4-5		± 1 kV Differential ± 2 kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Power Frequency 50/60 Hz Magnetic field IEC 61000-4-8		3 A/m	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

11.0 Notes on EMC

UK

Immunity Test	IEC 60601-Test Level	Compliance Level	Electromagnetic Environment - Guidance
Voltage Dips / Dropout IEC 61000-4-11		<p>< 5 % U_T (> 95 % Dip of the U_T) for 0.5 Cycle</p> <p>40 % U_T (60% Dip of the U_T) for 5 Cycles</p> <p>70% U_T (30 % Dip of the U_T) for 25 Cycles</p> <p>< 5 % U_T (>95 % Dip of the U_T) for 5 s</p>	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Varioair 3 demands continued function even in case of interruptions of the energy supply, it is recommended to supply the ATMOS Varioair 3 from an uninterruptible current supply or a battery.
NOTE U_T is the mains alternating current prior to application of the test levels.			

11.3 Guidelines and Manufacturer's Declaration - Immunity

The ATMOS Varioair 3 is intended for use in the electromagnetic environment specified below. The customer or user of the ATMOS Varioair 3 should ensure that it is used in such an environment.

Immunity Test	IEC 60601-Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	[3] V	Portable and mobile communications equipment should be separated from the ATMOS Varioair 3 incl. the cables by no less than the distances calculated/listed below. Recommended distances: $d = [1,2] \sqrt{P}$ $d = [1,2] \sqrt{P}$ $d = [2,3] \sqrt{P}$ where „P“ is the max. power in watts (W) and D is the recommended separation distance in meters (m). Field strengths from fixed transmitters, as determined by an electromagnetic site (a) survey, should be less than the compliance level (b). Interference may occur in the vicinity of equipment containing following symbol.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	[3] V/m	
			

11.0 Notes on EMC

UK

NOTE 1 With 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2 These guidelines don't like to be applicable in any case. The propagation of electromagnetic sizes is influenced by absorptions and reflections of buildings, objects and people.

a The field strength of stationary transmitters, such as base stations of cellular phones and mobile terrain radio equipment, amateur radio transmitters, cbm broadcast and TV stations cannot be predestined exactly. To determine the electromagnetic environment in regard to stationary transmitters, a study of the location is to be considered. If the measured field strength at the location where the Varioair 3 is used exceeds the above compliance level, the Varioair 3 is to be observed to verify the intended use. If abnormal performance characteristics are noted, additional measures might be necessary, e. g. a changed arrangement or another location for the device.

b Within the frequency range of 150 kHz to 80 MHz the field strength is to be below 3 V/m.

11.4 Recommended separations between portable and mobile RF Communications equipment and the ENT-camera Varioair 3

The Varioair 3 is intended for use in electromagnetic environment in which radiated disturbances are controlled. The customer or user of the Varioair 3 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications equipment and the Varioair 3 as recommended below, according to the maximum output power of the communications equipment.

Nominal output of the transmitter W	Separation distance, depending on transmit-frequency m		
	150 kHz bis 80 MHz $d = [1,2] \sqrt{P}$	80 MHz bis 800 MHz $d = [1,2] \sqrt{P}$	800 MHz bis 2,5 GHz $d = [2,3] \sqrt{P}$
0.1	0.12	0.12	0.23
0.1	0.46	0.46	0.9
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters for which the maximum nominal output isn't indicated in the above table, the recommended separation distance d in meters (m) can be determined using the equation belonging to the respective column whereas P is the maximum nominal output of the transmitter in watts (W) acc. to manufacturer's specification.

NOTE 1 With 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2 These guidelines don't like to be applicable in any case. The propagation of electromagnetic sizes is influenced by absorptions and reflections of buildings, objects and people.

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**EC - DECLARATION OF CONFORMITY
FOR MEDICAL PRODUCTS**

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Qd 148-4_CE0124

ATMOS General Standard Terms and Conditions

1. General:

Our General Standard Terms and Conditions apply exclusively. Client's terms and conditions which are contrary to or deviate from our General Standard Terms and Conditions are not recognised unless their validity is explicitly confirmed in writing. Our General Standard Terms and Conditions also apply even if we deliver to clients without reservation, in the knowledge of the client's contrary terms and conditions. Our General Standard Terms and Conditions also apply to all future business with that client.

2. Proposal - Order Confirmation

Our proposals are subject to change without notice unless otherwise stated in our order confirmation. Each order is only accepted by us following our written order confirmation.

3. Orders

Every order requires an exact description of all of our product's details. We assume no liability for errors and damage caused by inaccurate or incomplete ordering details.

4. Prices

Unless otherwise stated in the order confirmation, our prices in the order confirmation are ex factory prices and exclude packaging and value added tax. Packaging is charged separately at cost price in the invoice. Value added tax is charged separately in the invoice according to the legal rate on the invoice date. We reserve the right to change prices appropriately should price reductions or increases, especially due to wage settlements, changes in the price of materials or currency fluctuations, be incurred. Proof of such changes will be provided for the client on request.

5. Payment Conditions - Balancing

Unless otherwise stated in the order confirmation, our invoices are payable with a 3% discount within 10 days (except for repair and assembly services) or within 21 days from the invoice date net cash; money receipts is decisive for complying with this term. We are entitled to charge interest after the due date at a rate 2% above the relevant basic interest rate of the German Federal Bank. Should the client have payment arrears, we are entitled to charge interest on arrears at a rate 5% above the relevant basic interest rate of the German Federal Bank. Should we be able to prove higher damages due to arrears, we are also entitled to claim these. The client only has the right to balance invoices against its own claims should such claims be confirmed in a court of law or recognised by us. The client does not have the right of retention due to disputed counterclaims.

6. Delivery Periods

Fulfillment of our delivery duties requires the punctual and proper fulfilment of the client's duties. The right to defense on the grounds of an unfulfilled contract is reserved. Should the client default in accepting the goods delivery or breach other cooperation duties, we are entitled either to withdraw from the contract or claim compensation for any increased costs incurred up to that time without setting a further deadline. The right to make further claims is reserved. Furthermore, in such cases, the risk of coincidental destruction or a coincidental deterioration in the quality of the delivered goods is transferred to the client in the case of default in accepting such goods or payment arrears. Acts of God or stoppages (due to insufficient supplies of material, industrial disputes etc.) entitle us either to demand an appropriate extension of delivery periods or to partly or entirely dissolve the delivery contract. This does not give the client the right to claim damages. We have fulfilled delivery periods if the delivery goods have left our factory or the client has been informed of the goods' readiness for delivery within such delivery periods. Delivery periods stipulated by the client are not recognised by us unless they form part of our order confirmation. We adhere to legal terms and conditions in cases where, as a result of an undue delay in the delivery for which we are liable, the client is entitled to claim that his interests in a continued fulfilment of the contract have ceased. We also adhere to legal terms and conditions should a delay in delivery be caused by deliberate or grossly negligent action by us or our representatives for which we are responsible. We are also responsible for such actions by our representatives or agents. Should the delivery delay not be caused by our deliberate infringement of contractual duties for which we are responsible, our liability is limited to damage which is regarded as typical for that case. We are liable according to the legal terms and conditions if and in so far as the delivery delay for which we are responsible is caused by an infringement of a substantial contractual duty. In such cases, our liability is also limited to damage which is regarded as typical for that case. Should the delivery delay be caused by a culpable infringement of non-substantial contractual duties, our client is also entitled

to claim a one-off damage compensation worth 3 percentage points of the delivery value of the goods for each week's delay, up to a maximum which is no higher than 15 percentage points of the delivery value of the goods

7. Delivery - Familiarisation

In the case of the delivery of devices for the medico-technical industry which require assembly and/or familiarisation for the final customer using specialist trade personnel (such as Ear, Nose and Throat Apparatus and Suction Units), we reserve the right to deliver the goods exclusively to the relevant specialist traders. Should the trader not carry out assembly and/or familiarisation for the final customer, this is carried out by us. In such cases, we reserve the right to charge the client for the additionally created costs. Our specialist traders operate a recording system so that, if necessary, our products can be traced to the final customer. The specialist trader undertakes to immediately report to us all events and risks which must be reported in connection with our products.

8. Passage of Risk - Packaging

Unless otherwise stated in our order confirmation, delivery is agreed ex factory. The risk of the goods' damage or loss is therefore transferred to the client as soon as the goods leave the factory or the client is in default of acceptance of the goods. This also applies to cases where we confirm prepaid carriage. Transport packaging and all other packaging according to the packaging regulations is not returnable. Our client is responsible for disposing of the packaging at its own cost. Our deliveries are insured by us at the client's expense unless explicitly otherwise agreed. No insurance is arranged in the case of goods which are collected by our clients. In the case of transport damage, claims are only handled if the client receives confirmation of any damage, reduced weight or loss by the shipping company before accepting the delivery.

9. Warranty

The client is responsible for examining the delivered goods immediately after receiving them to determine any eventual deficiencies or delivery errors, and to report these immediately. Should the client fulfil this examining and reporting responsibility, and should payment conditions be fulfilled, we shall be liable to the client within the scope of legal regulations. Our period of warranty shall in all cases be two years. Our client can make use of the warranty as follows, so long as he can provide first buyer proof (in the form of an invoice or delivery note) and provided that the product still has the original, unchanged serial number:

- We choose whether to fulfil our guarantee by providing repair services free of charge - either on the client's premises or in our factory - or replacing the product. We can also provide these guarantee services through an authorised company;
- Should a product be returned to us, the client agrees to send the product in its original or similar packaging, offering the same protection as the original packaging, to our address or any address notified by us.
- Our guarantee ceases to apply if changes of any kind have been made to our product, unless such changes have been made by us or a company authorised by us, or have been previously agreed upon in writing by us. Our guarantee also ceases to apply if third parties have carried out repairs to our products or replaced parts thereof. This applies regardless of the fact whether these measures individually or collectively led to a deficiency of the product;
- We accept no responsibility for damage defects caused by
 - operational wear and tear;
 - incorrect installation or incorrect or insufficient maintenance;
 - incorrect operation of the product (in contradiction to the handbook delivered with the product);
 - improper use or operating faults;
 - inappropriate or negligent handling and care, especially with respect to dirt, lime, suction of fluids, inappropriate cleaning and sterilisation;
 - using accessories and/or replacement parts which are not explicitly approved;
 - incorrect assembly and/or initial operation by the client or third parties;
 - the client's negligence in handling the product;
 - unacceptable operating conditions, such as humidity, temperatures, the power supply, vibrations.
 - accidents, acts of God, especially lightning, water, fire, public unrest and insufficient ventilation. We are not liable for damage to other objects apart from our product itself, except in the case of any deliberate or grossly negligent actions by us or our representatives or agents. Should no deliberate breach of contract be claimed, our liability is limited to damage which is regarded as typical for that case. This also applies in the case of our culpable infringement of substantial contractual duties. The indispensable conditions of German Liability Law remain unaffected thereby.

- For second-hand equipment, the period of warranty shall be reduced to a period of twelve months.

10. Reservation of Ownership

We retain ownership of our goods until the receipt of all payments arising from the business relationship, including all demands arising from installation orders, subsequent orders, repairs, accessory deliveries and replacement orders. Should we have agreed upon payment on the basis of cheque and bill transactions, the ownership reservation applies until the cheque received by us has been paid in, and does not expire through our credit upon receiving the client's cheque. In the case of a breach of contract by the client, especially payment arrears, we are entitled to repossess our goods. Repossession of our goods represents a withdrawal from the contract, unless explicitly declared in writing by us. We have the right to utilise the product after its repossession, whilst the income form such use is balanced against the client's arrears, after deducting appropriate utilisation costs. The client is responsible for handling the goods with care. Should maintenance and inspection work be necessary, the client must carry these out punctually at his own cost. Our client is entitled to sell the goods he has bought from us in a proper sale transaction. However, he must immediately assign all outstanding claims to the value of the final invoice sum (including value added tax) of our claims to his customers or third parties. The client is entitled to collect this claim even after such assignment. Our right to collect the claim ourselves remains unaffected thereby. We undertake to release the securities to which we are entitled if requested to do so by the client should the realisable value of our securities be more than 10 percentage points higher than the outstanding claims. We reserve the right to choose the securities to be released.

11. Plans and Illustrations

We retain ownership of and copyrights to all plans, illustrations, calculations and other documents which are attached to our proposals. The client must receive explicit written permission before passing these on to third parties. Imitating our legally patented products is forbidden and will be prosecuted.

12. Jurisdiction and Place of Performance

Our central office is the place of performance for all disputes in connection with these General Standard Terms and Conditions and the contracts closed with clients under them. This jurisdiction excludes other jurisdiction relating to persons or subject-matter. Furthermore, our client is not entitled to bring charges against us in another court should he file counter-charges, carry out counterbalancing or declare retention. We, however, are entitled to bring charges against our client at their general place of jurisdiction or at another relevant court recognised by German or foreign law. Unless otherwise stated in the order confirmation, our central office is the place of performance.

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